

Mouse Brain Lysate

Cat. # BL7011

Size 100 μ l

Lysate Preparation:

Brain (cerebrum and cerebellum) was dissected from adult mouse (Balb/c). The tissue was rinsed in Hank's buffered saline solution and lysed in 1% SDS, 1.0 mM sodium ortho-vanadate, 10 mM Tris (pH 7.4) buffer. Brain lysates were sonicated then diluted to the final concentration in electrophoresis sample buffer.

Buffer and Storage:

Cell Lysates are supplied at a concentration of 1 mg/ml in electrophoresis sample buffer (62.5 mM Tris pH 6.8, 2% SDS, 5% glycerol, 0.003% bromophenol blue, 0.9% β -mercaptoethanol). Store at -20°C . Do not boil or dilute. Stable for 1 year.

Applications:

WB 20 μ l/lane

End user should determine optimal quantity for their particular applications and experiments.

Related Products:

LL7041 Mouse Liver Lysate
ML7051 Mouse Muscle Lysate
HL7021 Mouse Heart Lysate
SL7061 Mouse Spleen Lysate
KL7031 Mouse Kidney Lysate
MK6900 Mouse Tissue Lysate Sampler Kit

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

www.ecmbiosciences.com
telephone: 859-879-2075
toll-free: 1-800-859-8202
email: info@ecmbiosciences.com

ECMBiosciences

Rev 10/19/2007