

# $\alpha$ 1-Catenin (a.a. 143-153) Peptide

Cat. # CX3435

Size 50  $\mu$ g

## **Background:**

$\alpha$ -catenins are cadherin interacting proteins with homology to vinculin. Three  $\alpha$ -catenin genes have been described including  $\alpha$ 1-catenin ( $\alpha$ E-catenin),  $\alpha$ 2-catenin ( $\alpha$ N-catenin), and  $\alpha$ 3-catenin ( $\alpha$ T-catenin).  $\alpha$ 1-catenin has 81% homology with  $\alpha$ 2-catenin and 60% homology with  $\alpha$ 3-catenin. These  $\alpha$ -catenin isoforms may have similar roles since each binds cadherins. However, their expression patterns are both overlapping and distinct.  $\alpha$ 1-catenin was identified in epithelial cells, and is expressed in various cell types.  $\alpha$ 2-catenin is enriched in the nervous system, and  $\alpha$ 3-catenin is expressed highest in testis and heart. Phosphorylation may regulate the activity of  $\alpha$ 1-catenin, since tyrosine phosphorylation of Tyr-148 occurs during intercellular adhesion. This site is dephosphorylated by SHP2, which inhibits  $\alpha$ 1-catenin binding to  $\beta$ -catenin and translocation to the plasma membrane. Phosphorylation of  $\alpha$ 1-catenin at Tyr-148 may be important for inhibition of cell transformation, and dephosphorylation of this site may be important during SHP2-mediated cell transformation.

## **References**

- Herrenknecht, K. et al. (1991) Proc Natl Acad Sci U S A. 88(20):9156.  
Hirano, S. et al. (1992) Cell. 70(2):293.  
Janssens, B. et al. (2001) J Cell Sci. 114(17):3177.  
Burks, J. & Agazie, Y.M. (2006) Oncogene 25:7166.

## **Peptide Sequence:**

$\alpha$ 1-Catenin synthetic peptide corresponds to amino acid residues 143 to 153 in human  $\alpha$ 1-Catenin. This peptide sequence is highly conserved in rat and mouse  $\alpha$ 1-Catenin, and has some homology to  $\alpha$ 2-Catenin or  $\alpha$ 3-Catenin.

## **Buffer and Storage:**

Blocking Peptide is supplied in 50 $\mu$ l phosphate-buffered saline and 0.05% sodium azide. Store at  $-20^{\circ}$ C. Stable for 1 year.

## **Applications:**

Blocking 1:1000

ELISA 50 ng/well

End user should determine optimal concentration dependent on the concentration of the antibody.  
Recommended for blocking antibody reactivity in Western blot and immunocytochemistry.  
ELISA established in 96-well Nunc immunoplates where peptide was bound to plates for 2 hrs in 0.1 M sodium carbonate buffer, pH 8.5.

## **Specificity:**

The peptide is specifically recognized by  $\alpha$ 1-Catenin (a.a. 143-153) (CP3431) in ELISA, and has been shown to block the reactivity of CP3431 in Western blot. In addition, the peptide is recommended for use in blocking CP3431 reactivity in immunocytochemistry.

## **Related Products:**

- CP3431  $\alpha$ 1-Catenin (a.a. 143-153) Rabbit Polyclonal  
CP3451  $\alpha$ 1-Catenin (Tyr-148), phospho-specific Rabbit Polyclonal  
CX3455 phospho- $\alpha$ 1-Catenin (Tyr-148) Peptide  
CK6120  $\beta$ -Catenin Phospho-Regulation Antibody Sampler Kit  
CK6150  $\gamma$ -Catenin Phospho-Regulation Antibody Sampler Kit  
SM1631 SHP2 (N-terminal Region) Mouse Monoclonal

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

www.ecmbiosciences.com  
telephone: 859-879-2075  
toll-free: 1-800-859-8202  
email: info@ecmbiosciences.com

**ECM**Biosciences

Rev 6/29/2009